



# MEDISEPT

Perfect disinfection



## Velox® Spray

### Spray for surfaces' cleaning and disinfection

Alcohol liquid for rapid disinfection of small areas in the medical area, nursing homes, sanatoriums, fitness and SPA centers, hairdressers, massage and beauty salons, kindergartens and schools. The product is suitable for surfaces having and not-having contact with food.

#### Product features:

- For disinfection of small surfaces of medical equipment: beds and treatment armchairs, medical and surgical equipment, rehabilitation equipment.
- Designed for all non-medical surfaces e.g. in gastronomy, hairdressing and for surfaces having contact with food.
- Active already in 15 sec. against Fungi and in 30 sec. against Bacteria, Tuberculosis and viruses (including Rota, Noro).
- Does not leave streaks nor stains.
- Recommended especially for glass and shiny surfaces.

#### User benefits:

- Guaranteed effectiveness and safety of use.
- High reduction of the infection risk.
- Rapid effect of clean and disinfected surface.
- Can be used on many materials.



Non-invasive  
medical devices



Hairdressing  
and cosmetics  
tools



Tea tonic  
scent



Neutral  
scent

# Velox® Spray

## Application:

Alcohol liquid for rapid disinfection of small areas in the medical area, nursing homes, sanatoriums, fitness and SPA centers, hairdressers, massage and beauty salons, kindergartens and schools. The product is suitable for surfaces having and not-having contact with food. It can be used for disinfecting surfaces such as table tops, tables, handles, for hairdressing and cosmetics equipment and accessories. The product is recommended also for disinfection of non-invasive medical devices such as treatment chairs and medical equipment. Velox® Spray has bactericidal, yeasticidal activity against tuberculosis, limited virucidal activity (against Rotavirus, Norovirus) and against enveloped viruses (including Vaccinia, BVDV, HIV, HBV, HCV).

Effectiveness	EN 14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 13727	30 sec.	30 sec.
	EN 16615	30 sec.	30 sec.
Fungi ( <i>C. albicans</i> )	EN 13624	15 sec.	15 sec.
	EN 16615	30 sec.	30 sec.
<i>Mycobacteria tuberculosis</i> ( <i>M. terrae</i> )	EN 14348	30 sec.	30 sec.
Enveloped viruses ( <i>Vaccinia</i> , HIV, HBV, HCV, Herpes simplex, Ebola, SARS-Cov-2)	EN 14476	30 sec.	60 sec.
BVDV	EN 14476	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	60 sec.
Noro virus	EN 14476	30 sec.	30 sec.

## Directions for use:

Spray the surface from the distance of about 30 cm, ensuring it is fully covered. Leave it for 30 seconds. In case of disinfecting surfaces having contact with food, rinse them with drinking water. Other surfaces do not require rinsing with water. Do not use on surfaces sensitive to alcohol (e.g. acrylic glass). For professional use only.

## Ingredients: active substance in 100 g of the product:

63.7 g ethanol, 6.3 g propan-2-ol

## Packaging units:

500 ml spray bottle,  
1 L spray bottle,  
5 L canister



Medical device class II within the meaning of the Medical Devices Act. Intended for disinfecting non-invasive medical devices.  
Permission number for biocidal products in Poland: 7675/19

The safety data sheet and the component data sheet are available at [www.medisept.pl](http://www.medisept.pl)

CE 2274

 Producer:  
**MEDISEPT Sp. z o.o.**  
Konopnica 159 c  
21-030 Motycz

phone: + 48 81 535 22 22  
[export@medisept.pl](mailto:export@medisept.pl)  
[www.medisept.pl](http://www.medisept.pl)

Revision: 04.2020