



MEDISEPT

Perfect disinfection



Viruton® Forte

Concentrate for cleaning and disinfection of surgical instruments

Highly effective concentrate for simultaneous cleaning and disinfection of surgical and rotary instruments.

Product features:

- Concentration from 0.5% in 15 min.
- Active against Bacteria, Fungi, Tuberculosis, enveloped Viruses, Adeno and Polio.
- For tools made of steel, nickel, copper, aluminum.
- Cleans and flushes blood, protein and secretions.
- Anticorrosion properties and high material tolerance - Chirmed opinion (instruments manufacturer).
- Recommended for endoscopes - Vimex opinion (endoscopes manufacturer).
- For ultrasonic cleaners.

User benefits:

- Can be used on many sensitive materials.
- Guaranteed effectiveness and safety of use.
- High reduction of the infection risk.
- Effectively removes dried organic soiling.
- Economical.
- Speed and convenience in preparing the solution.



Concentrate



Active against B, F, V,
Tbc, Adeno, Polio



Instruments



Concentration
from 0.5%
in 15 min.

Viruton® Forte

Application:

Viruton® Forte has a very good cleaning properties and wide microbiological effectiveness within a short period of time. It has been designed for manual instrument disinfection and disinfection in ultrasonic cleaners. The product has anti-corrosion properties and a high material tolerance. Viruton® Forte is recommended for instruments made of stainless steel, nickel, copper, aluminum, rubber, porcelain, glass and plastics.

Effectiveness	EN 14885	Time	
		15 min.	30 min.
Bacteria (including MRSA)	EN 13727	0,5%*, 1%	
	EN 14561	1%	0,5%
	EN 13697	1%*	
Fungi (C. albicans)	EN 13624	0,5%	
	EN 14562	1%	0,5%
	EN 13697	1%*	
Mycobacteria (M. terrae, M. avium)	EN 14348, EN 14563	4%	1%
Enveloped viruses (Vaccinia, BVDV, HBV, HIV, HCV, SARS, Herpes Simplex, Influenza A, Ebola, SARS-Cov-2)	EN 14476	0,5%	
Adeno Virus	EN 14476	1%	
Polio Virus	EN 14476	4%	

* clean conditions

Directions for use:

Prepare working solution by diluting the concentrate with water. The level of dilution depends on a desired spectrum and duration of action. In order to obtain a concentration of 1% dilute 10 ml of the concentrate in 990 ml of water. Instruments shall be immersed in the solution immediately after use. All parts and holes must be completely filled in and covered by the solution. Once disinfected, rinse the instruments with water and dry. Activity of uncontaminated working solution is 14 days. For hygiene reasons, it is recommended to exchange the working solution on a daily basis. Do not combine with other detergents nor disinfectants. For professional use only.

Ingredients: 100 g contains:

3.76 g of N- (3-aminopropyl) -N-dodecylpropane-1,3-diamine (CAS: 2372-82-9), 3.39 g of Poly (oxy- 1,2-ethanedyl), alpha. [2 (didecylmethylamino) ethyl] - omega.-hydroxy-, propanoate (salt) (CAS: 94667-33-1), 0.75 g 2-aminoethanol (CAS: 141-43-5)

Packaging units:


1 L bottle with doser,
5 L canister



Medical device class II within the meaning of the Medical Devices Act. Intended for disinfection of non-invasive medical devices.

The safety data sheet and the component data sheet are available at www.medissept.pl

CE 2274

 Manufacturer:
MEDISEPT Sp. z o.o.
Konopnica 159 c
21-030 Motycz

phone: + 48 81 535 22 22
export@medissept.pl
www.medissept.pl

Revision: 04.2020