



# MEDISEPT

Perfect disinfection



## Velox® Wipes

### Alcohol wipes for cleaning and disinfection medical equipment's surfaces

Dedicated for rapid disinfection of small, difficult to access surfaces of medical equipment, chairs, beds, dental ends before the sterilisation process and prosthetic spoons and silicon moulds.

#### Product features:

- Active in 30 sec. against Bacteria, Fungi, Tuberculosis, enveloped viruses, Adeno and Rota virus.
- Recommended for units, medical equipment, rehabilitation equipment, dental handpieces, treatment lamps, patient beds.
- High material tolerance to ABS and upholstery materials - Famed opinion.
- Non-woven wipes 13x20 cm, 23 g / m<sup>2</sup>.

#### User benefits:

- Guaranteed effectiveness and safety of use.
- High reduction of infection risk.
- Rapid effect of clean and disinfected surface.
- Can be used on many sensitive materials.
- Refreshing effect.
- Economical to use.



Wipes



Active against  
Adeno, Rota



Non-invasive  
medical devices



Active  
in 30 sec.

## Application:

Presaturated with an alcohol-based formula, Velox® Wipes are a cleaner and disinfectant in one. Dedicated for rapid disinfection of small, difficult to access surfaces of medical equipment, chairs, beds, dental ends before the sterilisation process and prosthetic spoons and silicon moulds. Free from aldehydes and phenols, does not discolour disinfected surfaces. Thanks to an innovative formula based on the synergy of surfactants and nano-micelles, the product has a broad microbiological spectrum and unique cleaning properties. Tested and recommended by Famed - a leading European producer of top quality hospital equipment.

Effectiveness	EN14885	Conditions	
		clean	dirty
Bacteria (including MRSA)	EN 1276, EN 13697	30 sec.	60 sec.
	EN 13727	30 sec.	
	EN 16615	30 sec.	30 sec.
Fungi (C. albicans)	EN 1650	30 sec.	60 sec.
	EN 13624	30 sec.	
	EN 16615	30 sec.	30 sec.
Mycobacteria tuberculosis (M. terrae)	EN 14348	30 sec.	60 sec.
Enveloped viruses (HIV, HBV, HCV, influenza, Ebola, Vaccina, BVDV, SARS-Cov-2)	DVV/RKI	30 sec.	30 sec.
Rota virus	EN 14476	30 sec.	
Adeno virus	EN 14476	60 sec.	

## Directions for use:

Place the refill in the box. Take a tissue out of the box and wipe a surface to be cleaned and disinfected. Close the lid. After opening the package, use the tissues within 21 days. Do not use if a tissue is dry. Do not use on areas sensitive to alcohol (e.g. acrylic glass). For professional use only.

## Ingredients:

**100 pieces of tissues contain: 173 g of saturated formulation, 100 g contains:**

40 g pro-pan-2-ol (CAS: 67-63-0), 20 g ethanol, 0.39 g amine, n-C10-16 alkyltrimethylenedio di-reaction products with chloroacetic acid, 0.25 g didecyl-dimethylammonium chloride

## Packaging units:

100 pieces box,  
100 pieces refill



Medical device class II within the meaning of the Medical Devices Act. Dedicated for disinfecting non-invasive medical devices.

The safety data sheet and the component data sheet are available at [www.medisept.pl](http://www.medisept.pl)

**CE 2274**